




32222 Yamaha Tracer 900 Seat Rack Fitting Instructions

Kit Includes:

1X Seat Rack C32222

- The Pyramid seat rack simply fits in place of the passenger grabs rails.
- Remove the passenger seat and undo all the bolts shown in the diagram.
- Now refit the Pyramid seat rack to the same mounting points using the same fixings and refit the seat.
- Be sure to check fixings periodically.



Any questions?
sales@pyramid-plastics.co.uk
+44 (0)1427 677990
www.pyramid-plastics.co.uk

Follow us on social media and feel free to leave us a review!




32222 Yamaha Tracer 900 Seat Rack Fitting Instructions

Kit Includes:

1X Seat Rack C32222

- The Pyramid seat rack simply fits in place of the passenger grabs rails.
- Remove the passenger seat and undo all the bolts shown in the diagram.
- Now refit the Pyramid seat rack to the same mounting points using the same fixings and refit the seat.
- Be sure to check fixings periodically.



Any questions?
sales@pyramid-plastics.co.uk
+44 (0)1427 677990
www.pyramid-plastics.co.uk

Follow us on social media and feel free to leave us a review!

